



DIAM4 Server

DICOM Conformance Statement

Revision 2.2

07/04/2006

| Date | Rev. | Author | Comment |
|------------|------|---------------------------|---------------------------------------------------------|
| 10/02/2003 | 1 | H. Baekeland | Initial Release |
| 31/03/2005 | 2 | H. Baekeland | Added Query Retrieve SCU And Storage Commitment SCU/SCP |
| 28/04/2005 | 2.1 | H. Baekeland | Update Layout |
| 07/04/2006 | 2.2 | S. Moulinet JC. Judais | Added MPPS Added Presentation State |

Table of Contents

| | | |
|---------|----------------------------------------------------|----|
| 1 | Introduction | 3 |
| 1.1 | Purpose of this document | 3 |
| 1.2 | Sources for this document | 3 |
| 1.3 | Acronyms and abbreviation | 3 |
| 2 | Implementation Model | 4 |
| 2.1 | Application Data flow diagram | 5 |
| 2.2 | Functional Definition of Application Entities..... | 6 |
| 2.2.1 | Verification | 6 |
| 2.2.2 | Receive Images..... | 6 |
| 2.2.3 | Transmit Images | 6 |
| 2.2.4 | Query from Other Devices..... | 6 |
| 2.2.5 | Query to Other Devices | 6 |
| 2.2.6 | Retrieve to Other Devices | 6 |
| 2.2.7 | Retrieve from Other Devices | 6 |
| 2.2.8 | Receive Films..... | 7 |
| 2.2.9 | Send Films | 7 |
| 2.2.10 | Commitment to Store Images received..... | 7 |
| 2.2.11 | Ask Commitment to other devices..... | 7 |
| 2.2.12 | Modality Performed Procedure Step | 7 |
| 2.2.13 | Presentation state | 7 |
| 2.3 | Sequencing of real world activity | 7 |
| 3 | AE Specifications..... | 8 |
| 3.1 | DIAM ⁴ Server specification | 8 |
| 3.1.1 | Association establishment Policies | 11 |
| 3.1.1.1 | General..... | 11 |
| 3.1.1.2 | Number of Associations | 11 |
| 3.1.1.3 | Asynchronous nature | 11 |
| 3.1.1.4 | Implementation Identifying Information..... | 11 |
| 3.1.2 | Association Initiation by Real World Activity..... | 12 |
| 3.1.2.1 | Real World Activity – Find | 12 |
| 3.1.2.2 | Real World Activity - Move Images..... | 13 |
| 3.1.2.3 | Real World Activity - Storage as SCU..... | 14 |
| 3.1.2.4 | Real World Activity – Storage Commitment | 15 |
| 3.1.2.5 | Real World Activity – Print as SCU..... | 17 |
| 3.1.2.6 | Real World Activity – MPPS as SCU..... | 18 |
| 3.1.3 | Association Acceptance Policy..... | 18 |
| 3.1.3.1 | Real World Activity – Verification | 18 |
| 3.1.3.2 | Real World Activity - Storage as SCP | 19 |
| 3.1.3.3 | Real World Activity – Find | 21 |
| 3.1.3.4 | Real World Activity – Move..... | 25 |
| 3.1.3.5 | Real World Activity – Storage Commitment | 27 |
| 3.1.3.6 | Real World Activity – MPPS as SCP | 29 |
| 3.1.4 | Presentation State Display | 33 |
| 4 | Communications Profiles | 38 |
| 4.1 | TCP/IP Stack..... | 38 |
| 4.1.1 | Physical Media Support..... | 38 |
| 5 | Extensions/Specializations/Privatizations | 38 |
| 6 | Configuration | 38 |
| 7 | Support for Extended Character Sets | 38 |

1 Introduction

1.1 Purpose of this document

The purpose of this document is to describe how DIAM⁴ Server conforms to the DICOM standard. It describes what parts and definition it uses and in what way, in order to provide interoperability with other devices that claim same conformance.

1.2 Sources for this document

American College of Cardiology – National Manufactures Association (ACR-NEMA) Digital Imaging and Communications V2.0
ACR-NEMA Digital Imaging and Communications in Medicine (DICOM) v3.0, Final Draft, May. 1998

1.3 Acronyms and abbreviation

The following symbols and abbreviations are used in this Part.


| | |
|-------------|----------------------------------------------------------------------------------|
| ➤ ACR | American College of Radiology |
| ➤ ACSE | Association Control Service Element |
| ➤ AE | Application Entity |
| ➤ ANSI | American National Standards Institute |
| ➤ AP | Application Profile |
| ➤ API | Application Programming Interface |
| ➤ ASCII | American Standard Code for Information Interchange |
| ➤ CEN TC251 | Comite Europeen de Normalisation - Technical Committee 251 – Medical Informatics |
| ➤ DICOM | Digital Imaging and Communications in Medicine |
| ➤ DIMSE | DICOM Message Service Element |
| ➤ DIMSE-C | DICOM Message Service Element-Composite |
| ➤ DIMSE-N | DICOM Message Service Element-Normalized |
| ➤ FSC | File-set Creator |
| ➤ FSR | File-set Reader |
| ➤ FSU | File-set Updater |
| ➤ HISPP | Healthcare Informatics Standards Planning Panel |
| ➤ HL7 | Health Level 7 |
| ➤ IE | Information Entity |
| ➤ IEEE | Institute of Electrical and Electronics Engineers |
| ➤ IOD | Information Object Definition |
| ➤ ISO | International Standards Organization |
| ➤ ISP | International Standardized Profile |
| ➤ JIRA | Japanese Industry Radiology Apparatus |
| ➤ MSDS | Healthcare Message Standard Developers Sub-Committee |
| ➤ NEMA | National Electrical Manufacturers Association |
| ➤ OSI | Open Systems Interconnection |
| ➤ PDU | Protocol Data Unit |
| ➤ RWA | Real-World Activity |
| ➤ SCP | Service Class Provider |
| ➤ SCU | Service Class User |
| ➤ SOP | Service-Object Pair |

- TCP/IP Transmission Control Protocol/Internet Protocol
- UID Unique Identifier

2 Implementation Model

DIAM⁴ Server is a web based PACS solution. DIAM⁴ Server stores images sent to it by service class users, takes responsibility for storage of the images, can respond to a storage commitment demand, allows queries based on several standard query models, and retrieves requested images. DIAM⁴ Server also stores films sent to it by service class user.

2.1 Application Data flow diagram



2.2 Functional Definition of Application Entities

2.2.1 Verification

DIAM⁴ Server will respond to a **C-ECHO** verification.
DIAM⁴ Server will also ask for a **C-ECHO** verification.

2.2.2 Receive Images

DIAM⁴ Server stores a received image in its entirety, without change, in its internal data store. DIAM⁴ Server stores each image with the File Meta Information attached to it.

DIAM⁴ Server extracts the query information with respect to the patient, study, series and image, and stores this information within its internal database.

2.2.3 Transmit Images

DIAM⁴ Server acts a Service Class User of C-Store to transmit images to other compatible devices.

2.2.4 Query from Other Devices

DIAM⁴ Server responds to queries based on the records stored in its database.

2.2.5 Query to Other Devices

DIAM⁴ Server can query other compatible devices.

2.2.6 Retrieve to Other Devices

DIAM⁴ Server acts as a Service Class Provider of C-Move to retrieve images. It does so by obtaining a reference from the database then obtaining the image object itself from the data store.

2.2.7 Retrieve from Other Devices

DIAM⁴ Server can retrieve images from other compatible devices.

2.2.8 Receive Films

DIAM⁴ Server stores a received film in its entirety, without change, in its internal data store. DIAM⁴ Server stores each film with the File Meta Information attached to it.

2.2.9 Send Films

DIAM⁴ Server can send films to other compatible devices.

2.2.10 Commitment to Store Images received

DIAM⁴ Server acts a Service Class Provider of Storage Commitment to explicitly take responsibility for storing images received.

2.2.11 Ask Commitment to other devices

DIAM⁴ Server can ask storage commitment to other compatible devices.

2.2.12 Modality Performed Procedure Step

DIAM⁴ Server can receive Modality Performed Procedure Step messages from modalities. When DIAM⁴ Server receives a DICOM Modality Performed Procedure Step SOP Class N-CREATE request, it creates a new object in its database. When DIAM⁴ Server receives a DICOM Modality Performed Procedure Step SOP Class N-SET request, it updates the corresponding object.

DIAM⁴ Server can also route this messages to other Dicom devices.

2.2.13 Presentation state

DIAM⁴ Server can receive Grayscale Softcopy Presentation State objects from other DICOM devices. It stores these in its database and applies them when the display of related objects is performed.

DIAM⁴ Server can't create Grayscale Softcopy Presentation State objects.

2.3 Sequencing of real world activity

Not applicable.

3 AE Specifications

3.1 DIAM⁴ Server specification

DIAM⁴ Server provides Standard Conformance to the following DICOM V3.0 SOP Class as an SCU.

Table 1 : Verification SOP Class as SCU

| SOP Class Name | SOP Class UID |
|----------------|-------------------|
| Verification | 1.2.840.10008.1.1 |

Table 2 : Query/Retrieve SOP Classes as SCU

| SOP Class Name | SOP Class UID |
|-----------------------------------|-----------------------------|
| Study Root Query/Retrieve IM Find | 1.2.840.10008.5.1.4.1.2.2.1 |
| Study Root Query/Retrieve IM Move | 1.2.840.10008.5.1.4.1.2.2.2 |

Table 3 : Storage SOP Classes as SCU

| SOP Class Name | SOP Class UID |
|---------------------------------------------------------|-------------------------------|
| Basic Text SR Storage | 1.2.840.10008.5.1.4.1.1.88.11 |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 |
| Comprehensive SR Storage | 1.2.840.10008.5.1.4.1.1.88.33 |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 |
| Digital Intra Oral X-ray Image Storage For Presentation | 1.2.840.10008.5.1.4.1.1.1.3 |
| Digital Intra Oral X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.3.1 |
| Digital Mammography X-ray Image Storage For Present. | 1.2.840.10008.5.1.4.1.1.1.2 |
| Digital Mammography X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.2.1 |
| Digital X-ray Image Storage For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 |
| Digital X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 |
| Enhanced MR Image Storage | 1.2.840.10008.5.1.4.1.1.4.1 |
| Enhanced SR Storage | 1.2.840.10008.5.1.4.1.1.88.22 |
| Grayscale Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.1 |
| Hardcopy Color Image Storage | 1.2.840.10008.5.1.1.30 |
| Hardcopy Grayscale Image Storage | 1.2.840.10008.5.1.1.29 |
| Key Object Selection Document | 1.2.840.10008.5.1.4.1.1.88.59 |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 |
| MR Spectroscopy Storage | 1.2.840.10008.5.1.4.1.1.4.2 |
| PET Image Storage | 1.2.840.10008.5.1.4.1.1.128 |
| Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.20 |
| RETIRED Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.5 |
| RETIRED Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6 |
| RETIRED Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3 |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 |
| Stored Print Storage | 1.2.840.10008.5.1.1.27 |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 |
| Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 |
| X-ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 |
| X-ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 |

Table 4 : Storage Commitment SOP Classes as SCU

| SOP Class Name | SOP Class UID |
|-------------------------------|----------------------|
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 |

Table 5 : Print SOP Class as SCU

| SOP Class Name | SOP Class UID |
|------------------------|-------------------------|
| BasicFilmSession | 1.2.840.10008.5.1.1.1 |
| BasicFilmBox | 1.2.840.10008.5.1.1.2 |
| BasicGrayscaleImageBox | 1.2.840.10008.5.1.1.4 |
| BasicColorImageBox | 1.2.840.10008.5.1.1.4.1 |
| Printer | 1.2.840.10008.5.1.1.16 |

DIAM⁴ Server provides Standard Conformance to the following DICOM V3.0 SOP Class as an SCP.

Table 6 : Verification SOP Class as SCP

| SOP Class Name | SOP Class UID |
|----------------|-------------------|
| Verification | 1.2.840.10008.1.1 |

Table 7 : Query/Retrieve SOP Classes as SCP

| SOP Class Name | SOP Class UID |
|-------------------------------------------|-----------------------------|
| Patient Root Query/Retrieve IM Find | 1.2.840.10008.5.1.4.1.2.1.1 |
| Patient Root Query/Retrieve IM Move | 1.2.840.10008.5.1.4.1.2.1.2 |
| Study Root Query/Retrieve IM Find | 1.2.840.10008.5.1.4.1.2.2.1 |
| Study Root Query/Retrieve IM Move | 1.2.840.10008.5.1.4.1.2.2.2 |
| Patient/Study Root Query/Retrieve IM Find | 1.2.840.10008.5.1.4.1.2.3.1 |
| Patient/Study Root Query/Retrieve IM Move | 1.2.840.10008.5.1.4.1.2.3.2 |

Table 8 : Storage Commitment SOP Classes as SCP

| SOP Class Name | SOP Class UID |
|-------------------------------|----------------------|
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 |

Table 9 : Storage SOP Classes as SCP

| SOP Class Name | SOP Class UID |
|---------------------------------------------------------|-------------------------------|
| Basic Text SR Storage | 1.2.840.10008.5.1.4.1.1.88.11 |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 |
| Comprehensive SR Storage | 1.2.840.10008.5.1.4.1.1.88.33 |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 |
| Digital Intra Oral X-ray Image Storage For Presentation | 1.2.840.10008.5.1.4.1.1.1.3 |
| Digital Intra Oral X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.3.1 |
| Digital Mammography X-ray Image Storage For Present. | 1.2.840.10008.5.1.4.1.1.1.2 |
| Digital Mammography X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.2.1 |
| Digital X-ray Image Storage For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 |
| Digital X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 |
| Enhanced MR Image Storage | 1.2.840.10008.5.1.4.1.1.4.1 |
| Enhanced SR Storage | 1.2.840.10008.5.1.4.1.1.88.22 |
| Grayscale Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.1 |
| Hardcopy Color Image Storage | 1.2.840.10008.5.1.1.30 |
| Hardcopy Grayscale Image Storage | 1.2.840.10008.5.1.1.29 |
| Key Object Selection Document | 1.2.840.10008.5.1.4.1.1.88.59 |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 |
| MR Spectroscopy Storage | 1.2.840.10008.5.1.4.1.1.4.2 |
| PET Image Storage | 1.2.840.10008.5.1.4.1.1.128 |
| Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.20 |
| RETIRED Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.5 |
| RETIRED Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6 |
| RETIRED Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3 |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 |
| Stored Print Storage | 1.2.840.10008.5.1.1.27 |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 |
| Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 |
| X-ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 |
| X-ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 |

Table 10 : Print SOP Class as SCP

| SOP Class Name | SOP Class UID |
|------------------------|-------------------------|
| BasicFilmSession | 1.2.840.10008.5.1.1.1 |
| BasicFilmBox | 1.2.840.10008.5.1.1.2 |
| BasicGrayscaleImageBox | 1.2.840.10008.5.1.1.4 |
| BasicColorImageBox | 1.2.840.10008.5.1.1.4.1 |
| Printer | 1.2.840.10008.5.1.1.16 |

3.1.1 Association establishment Policies

3.1.1.1 General

The following Application Context Name will be proposed and recognized by DIAM4 Server

- DICOM 3.0 Application Context 1.2.840.10008.3.1.1.1

3.1.1.2 Number of Associations

The maximum number of association accepted or maintained by DIAM4 Server is limited only by the physical memory of the machine on which it runs. Typically it can be up to 20.

3.1.1.3 Asynchronous nature

DIAM⁴ Server allows a single outstanding operation on any association. Therefore, DIAM⁴ Server does not support asynchronous operations window negotiation, other than the default as specified by the DICOM specification.

3.1.1.4 Implementation Identifying Information

DIAM⁴ Server will respond with the following implementation identifying parameters:

Implementation Class UID 1.2.826.0.1.3680043.2.406.0001

The last number of the implementation class UID is the machine serial number.

Implementation Version Name GIOL_DIAM_400

The implementation version name policies are the following: product name "GIOL_DIAM_" followed by the version of the product, "400" (meaning "version 4.0.0").

3.1.2 Association Initiation by Real World Activity

3.1.2.1 Real World Activity – Find

3.1.2.1.1 Associated Real World Activity – Find

DIAM⁴ Server will issue a FIND request when a user of DIAM⁴ Server wishes to view patient and study information from a remote DICOM SCP.

3.1.2.1.2 Presentation context Table – Find

DIAM⁴ Server supports the transfer syntaxes listed in Table 11. For a QUERY request, DIAM⁴ Server supports the Presentation Contexts listed in Table 12.

Table 11 : Query Transfer Syntaxes

| Transfer Syntax | UID |
|-------------------------------------------------|---------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |
| DICOM Explicit VR Big Endian Transfer Syntax | 1.2.840.10008.1.2.2 |

Table 12 : Query Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|--------------------------------------|---------------------------------|----------------------|------|----------------------|
| SOP Class | SOP Class UID | | | |
| Study Root Query/Retrieve IM Find | 1.2.840.10008.5.1.4. 1.2.2.1 | all from Table 11 | SCU | None |

3.1.2.1.3 SOP Specific conformance – Find

DIAM⁴ Server supports the following search keys.

Table 13 : Patient level attributes

| Description | Tag |
|--------------|------------------|
| Patient name | (0x0010, 0x0010) |
| Patient id | (0x0010, 0x0020) |

Table 14 : Study level attributes

| Description | Tag |
|----------------------------|------------------|
| Study id | (0x0020, 0x0010) |
| Study date | (0x0008, 0x0020) |
| Study time | (0x0008, 0x0010) |
| Accession number | (0x0008, 0x0050) |
| Modalities in Study | (0x0008, 0x0061) |
| Referring Physician's Name | (0x0008, 0x0090) |
| Study Description | (0x0008, 0x1030) |
| Study Instance UID | (0x0020, 0x000D) |

Table 15 : Series level attributes

| Description | Tag |
|-------------------------------------|------------------|
| Series instance UID | (0x0020, 0x000E) |
| Series number | (0x0020, 0x0011) |
| Modality | (0x0008, 0x0060) |
| Series Description | (0x0008, 0x103e) |
| Request Attribute Sequence | (0x0040, 0x0275) |
| > Requested Procedure ID | (0x0040, 0x1001) |
| > Scheduled Procedure Step ID | (0x0040, 0x0009) |
| Performed Procedure Step Start Date | (0x0040, 0x0244) |
| Performed Procedure Step Start Time | (0x0040, 0x0245) |

Table 16 : Composite Object Instance Level attributes

| Description | Tag |
|------------------|------------------|
| SOP instance UID | (0x0008, 0x0018) |
| Instance Number | (0x0020, 0x0013) |

Table 17 : SR Document Specific Level attributes

| Description | Tag |
|-----------------------------|------------------|
| Completion Flag | (0x0040, 0xA491) |
| Verification Flag | (0x0040, 0xA493) |
| Verifying Observer Sequence | (0x0040, 0xA073) |
| > Verification DateTime | (0x0040, 0xA030) |
| > Verifying Observer Name | (0x0040, 0xA075) |
| Concept Name Code Sequence | (0x0040, 0xA043) |
| > Code Value | (0x0008, 0x0100) |
| > Coding Scheme Designator | (0x0008, 0x0102) |

Table 18 : Key Object Document Specific Level attributes

| Description | Tag |
|----------------------------|------------------|
| Concept Name Code Sequence | (0x0040, 0xA043) |
| > Code Value | (0x0008, 0x0100) |
| > Coding Scheme Designator | (0x0008, 0x0102) |

Table 19 : Presentation State Specific Level attributes

| Description | Tag |
|-------------------------------|------------------|
| Content Label | (0x0070, 0x0080) |
| Content Description | (0x0070, 0x0081) |
| Presentation Creation Date | (0x0070, 0x0082) |
| Presentation Creation time | (0x0070, 0x0083) |
| Presentation Creator's Name | (0x0070, 0x0084) |
| Referenced Series Sequence | (0x0008, 0x1115) |
| > Series Instance UID | (0x0020, 0x000E) |
| > Referenced Image Sequence | (0x0008, 0x1140) |
| >>Referenced SOP Class UID | (0x0008, 0x1150) |
| >>Referenced SOP Instance UID | (0x0008, 0x1155) |

3.1.2.2 Real World Activity - Move Images

3.1.2.2.1 Associated Real World Activity – Move Images

DIAM⁴ Server will issue a MOVE request when a user of DIAM⁴ Server wishes to move one or more studies / series from a remote DICOM SCP back to DIAM⁴ Server (retrieve) or another remote DICOM SCP.

3.1.2.2.2 Presentation context Table – Move

DIAM⁴ Server supports the transfer syntaxes listed in Table 20. For a MOVE request, DIAM⁴ Server supports the Presentation Contexts listed in Table 21.

Table 20 : Move Transfer Syntaxes

| Transfer Syntax | UID |
|-------------------------------------------------|---------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |
| DICOM Explicit VR Big Endian Transfer Syntax | 1.2.840.10008.1.2.2 |

Table 21 Move Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|-----------------------------------|-----------------------------|------------------------|-------------|-----------------------------|
| SOP Class | SOP Class UID | | | |
| Study Root Query/Retrieve IM Move | 1.2.840.10008.5.1.4.1.2.2.2 | all from Table 20 | SCU | None |

3.1.2.2.3 SOP Specific Conformance – Move

See SOP Specific Conformance – Storage SCU.

3.1.2.3 Real World Activity - Storage as SCU

3.1.2.3.1 Associated Real World Activity – Storage as SCU

DIAM⁴ Server will issue a Storage request when a user of DIAM⁴ wishes to send a study of images to a remote DICOM SCP.

3.1.2.3.2 Presentation context Table – Storage as SCU

DIAM⁴ Server supports the transfer syntaxes listed in Table 22. For a Storage request, DIAM⁴ Server supports the Presentation Contexts listed in Table 23.

Table 22 : Storage Transfer Syntaxes

| Transfer Syntax | UID |
|------------------------------------------------------------------------|------------------------|
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Big Endian Transfer Syntax | 1.2.840.10008.1.2.2 |
| DICOM Lossy JPEG 8 Bit – JPEG Baseline (Process 1) | 1.2.840.10008.1.2.4.50 |
| DICOM JPEG Lossless, hierarchical, first order prediction (Process 14) | 1.2.840.10008.1.2.4.70 |
| RLE Lossless Transfer Syntax | 1.2.840.10008.1.2.5 |

3.1.2.3.3 SOP Specific Conformance – Storage as SCU

Table 23 : Storage Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|------------------------|-------------------------------|------------------------|-------------|-----------------------------|
| SOP Class | SOP Class UID | | | |
| Basic Text SR Storage | 1.2.840.10008.5.1.4.1.1.88.11 | all from Table20 | SCU | None |

| | | | | |
|---------------------------------------------------------|--------------------------------|------------------|-----|------|
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | all from Table20 | SCU | None |
| Comprehensive SR Storage | 1.2.840.10008.5.1.4.1.1.88.33 | all from Table20 | SCU | None |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | all from Table20 | SCU | None |
| Digital Intra Oral X-ray Image Storage For Presentation | 1.2.840.10008.5.1.4.1.1.1.3 | all from Table20 | SCU | None |
| Digital Intra Oral X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.3.1 | all from Table20 | SCU | None |
| Digital Mammography X-ray Image Storage For Present. | 1.2.840.10008.5.1.4.1.1.1.2 | all from Table20 | SCU | None |
| Digital Mammography X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.2.1 | all from Table20 | SCU | None |
| Digital X-ray Image Storage For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | all from Table20 | SCU | None |
| Digital X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 | all from Table20 | SCU | None |
| Enhanced MR Image Storage | 1.2.840.10008.5.1.4.1.1.4.1 | all from Table20 | SCU | None |
| Enhanced SR Storage | 1.2.840.10008.5.1.4.1.1.88.22 | all from Table20 | SCU | None |
| Grayscale Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.1 | all from Table20 | SCU | None |
| Hardcopy Color Image Storage | 1.2.840.10008.5.1.1.30 | all from Table20 | SCU | None |
| Hardcopy Grayscale Image Storage | 1.2.840.10008.5.1.1.29 | all from Table20 | SCU | None |
| Key Object Selection Document | 1.2.840.10008.5.1.4.1.1.88.59 | all from Table20 | SCU | None |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 | all from Table20 | SCU | None |
| MR Spectroscopy Storage | 1.2.840.10008.5.1.4.1.1.4.2 | all from Table20 | SCU | None |
| PET Image Storage | 1.2.840.10008.5.1.4.1.1.128 | all from Table20 | SCU | None |
| Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.20 | all from Table20 | SCU | None |
| RETIRED Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.5 | all from Table20 | SCU | None |
| RETIRED Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6 | all from Table20 | SCU | None |
| RETIRED Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3 | all from Table20 | SCU | None |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | all from Table20 | SCU | None |
| Stored Print Storage | 1.2.840.10008.5.1.1.27 | all from Table20 | SCU | None |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | all from Table20 | SCU | None |
| Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | all from Table20 | SCU | None |
| X-ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | all from Table20 | SCU | None |
| X-ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 | all from Table20 | SCU | None |
| Raw data | 1.2.840.10008.5.1.4.1.1.66 | all from Table20 | SCU | None |
| VL photographic | 1.2.840.10008.5.1.4.1.1.77.1.4 | all from Table20 | SCU | None |

3.1.2.4 Real World Activity – Storage Commitment

3.1.2.4.1 Associated Real World Activity – Storage Commitment

DIAM⁴ Server stores images that are sent to it from an SCU. In some configurations (eg cache-only), DIAM⁴ Server may send images to another SCP (eg a PACS) for permanent storage. The request for storage commitment may then be transmitted from DIAM⁴ Server together with a list of references to one or more SOP instances. This action is invoked through the DIMSE N-ACTION primitive. The following message is supported:

- Request Storage Commitment - to request the safekeeping of a set of SOP instances

3.1.2.4.2 Presentation Context Table – Storage Commitment

DIAM⁴ Server supports the transfer syntaxes listed in Table 24. For a Storage Commitment request, DIAM⁴ Server supports the Presentation Contexts listed in Table 25.

Table 24 : Storage Commitment Transfer Syntaxes

| Transfer Syntax | UID |
|-------------------------------------------------|---------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |
| DICOM Explicit VR Big Endian Transfer Syntax | 1.2.840.10008.1.2.2 |

Table 25 : Storage Commitment Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|---------------------------------------------------|----------------------------------------------|------------------------|-------------|-----------------------------|
| SOP Class Storage Commitment Push Model | SOP Class UID 1.2.840.10008.1.20.1 | all from Table 23 | SCU | None |

3.1.2.4.3 SOP Specific Conformance – Storage Commitment

DIAM⁴ Server provides standard conformance to the DICOM **Storage Commitment** Service Class.

DIAM⁴ Server supports the following elements for this SOP class as a SCU. The Transaction UID Attribute (0008,1195) value generated by DIAM⁴ Server uniquely identifies each Storage Commitment Request.

Table 26 :Storage Commitment Request – Action Information

| Action Type Name | Action Type ID | Attribute Name | Tag |
|----------------------------|----------------|------------------------------|-------------|
| Request Storage Commitment | 1 | Transaction UID | (0008,1195) |
| | | Referenced SOP Sequence | (0008,1199) |
| | | >Referenced SOP Class UID | (0008,1150) |
| | | >Referenced SOP Instance UID | (0008,1155) |

Subsequently, DIAM⁴ Server expects N-EVENT-REPORT from the SCP. DIAM⁴ Server returns an N-EVENT-REPORT response primitive.

3.1.2.5 Real World Activity – Print as SCU

3.1.2.5.1 Associated Real World Activity – Print as SCU

DIAM⁴ Server will issue a Print request when a user of DIAM⁴ wishes to send images to a remote DICOM Printer SCP.

3.1.2.5.2 Presentation context Table – Print as SCU

DIAM⁴ Server supports the transfer syntaxes listed in Table 27 for a Print request; DIAM⁴ Server supports the Presentation contexts listed in Table 27.

Table 27 : Print Transfer Syntaxes

| Transfer Syntax | UID |
|-------------------------------------------------|---------------------|
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |

3.1.2.5.3 SOP specific Conformance – Print as SCU

Table 28 : Print Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|------------------------|----------------------------------|-------------------|------|----------------------|
| SOP Class | SOP Class UID | | | |
| BasicFilmSession | 1.2.840.10008.5.1.1.1 | all from Table 27 | SCU | None |
| BasicFilmBox | 1.2.840.10008.5.1.1.2 | all from Table 27 | SCU | None |
| BasicGrayscaleImageBox | 1.2.840.10008.5.1.1.4 | all from Table 27 | SCU | None |
| BasicColorImageBox | 1.2.840.10008.5.1.1.4 | all from Table 27 | SCU | None |
| Printer | .1 1.2.840.10008.5.1.1.1 6 | all from Table 27 | SCU | None |

3.1.2.6 Real World Activity – MPPS as SCU

3.1.2.6.1 Associated Real World Activity – MPPS SCU

DIAM⁴ Server will issue a Modality Performed Procedure Step request each time it receives valid MPPS N-CREATE or N-SET messages from other DICOM systems. These messages are routed to all MPPS destinations configured in DIAM⁴ Server.

3.1.2.6.2 Presentation context Table – MPPS SCU

DIAM⁴ Server will accept any of the transfer syntaxes listed in Table 55. DIAM⁴ Server will accept any of the Presentation Contexts listed in Table 56.

3.1.2.6.3 SOP Specific Conformance – MPPS SCU

DIAM⁴ Server supports the elements MPPS N-CREATE request as an SCU, contained in table 57 ; it also supports the elements MPPS N-SET request as an SCU, contained in table 58.

3.1.3 Association Acceptance Policy

3.1.3.1 Real World Activity – Verification

3.1.3.1.1 Associated Real World Activity – Verification

DIAM⁴ Server will respond to Verification requests provided by SCU with the ability to determine if DIAM⁴ Server can receive DICOM requests.

3.1.3.1.2 Presentation Context Table – Verification

DIAM⁴ Server supports the transfer syntaxes listed in Table 29. DIAM⁴ Server will accept any of the Presentation Contexts listed in Table 30 for Verification.

Table 29 : Verification Transfer Syntaxes

| Transfer Syntax | UID |
|-------------------------------------------------|---------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |
| DICOM Explicit VR Big Endian Transfer Syntax | 1.2.840.10008.1.2.2 |

Table 30 : Verification Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|-----------------|-------------------|-------------------|------|----------------------|
| SOP Class | SOP Class UID | | | |
| Verification | 1.2.840.10008.1.1 | all from Table 29 | SCP | None |

3.1.3.1.3 SOP Specific Conformance – Verification

DIAM⁴ Server provides standard conformance to the DICOM Verification Service Class. DIAM⁴ Server returns one of the following status codes.

Table 31 : Verification status codes.

| Service Status | Further Meaning | Protocol Codes | Related Fields | Description |
|----------------|-----------------|----------------|----------------|-------------|
|----------------|-----------------|----------------|----------------|-------------|

| | | | | |
|---------|---------|------|--|-----------------------------------|
| Error | Failed | C000 | | The operation was not successful. |
| Success | Success | 0000 | | Operation performed properly. |

3.1.3.1.4 Presentation Context Acceptance Criterion – Verification
 DIAM⁴ Server will always accept a Presentation Context for the Verification SOP Class with the default DICOM transfer syntax listed in Table 29.

3.1.3.1.5 Transfer Syntax Selection Policies – Verification
 Since no DICOM data object is associated with a Verification command, only the default DICOM transfer syntax is required/supported.

3.1.3.2 Real World Activity - Storage as SCP

3.1.3.2.1 Associated Real World Activity – Storage as SCP
 DIAM⁴ Server will archive images that are sent to it from an SCU.

3.1.3.2.2 Presentation Context Table – Storage as SCP
 DIAM4 Server supports the following transfer syntaxes listed in Table 32. DIAM4 Server supports any of the Presentation Contexts listed in Table 33 for Storage.

Table 32 : Storage Transfer Syntaxes

| Transfer Syntax | UID |
|------------------------------------------------------------------------|------------------------|
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Big Endian Transfer Syntax | 1.2.840.10008.1.2.2 |
| DICOM Lossy JPEG 8 Bit – JPEG Baseline (Process 1) | 1.2.840.10008.1.2.4.50 |
| DICOM JPEG Lossless, hierarchical, first order prediction (Process 14) | 1.2.840.10008.1.2.4.70 |
| RLE Lossless Transfer Syntax | 1.2.840.10008.1.2.5 |

Table 33 : Storage Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|---------------------------------------------------------|-------------------------------|-------------------|------|----------------------|
| SOP Class | SOP Class UID | | | |
| Basic Text SR Storage | 1.2.840.10008.5.1.4.1.1.88.11 | all from Table 23 | SCP | None |
| CT Image Storage | 1.2.840.10008.5.1.4.1.1.2 | all from Table 23 | SCP | None |
| Comprehensive SR Storage | 1.2.840.10008.5.1.4.1.1.88.33 | all from Table 23 | SCP | None |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | all from Table 23 | SCP | None |
| Digital Intra Oral X-ray Image Storage For Presentation | 1.2.840.10008.5.1.4.1.1.1.3 | all from Table 23 | SCP | None |
| Digital Intra Oral X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.3.1 | all from Table 23 | SCP | None |
| Digital Mammography X-ray Image Storage For Present. | 1.2.840.10008.5.1.4.1.1.1.2 | all from Table 23 | SCP | None |
| Digital Mammography X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.2.1 | all from Table 23 | SCP | None |
| Digital X-ray Image Storage For Presentation | 1.2.840.10008.5.1.4.1.1.1.1 | all from Table 23 | SCP | None |
| Digital X-ray Image Storage For Processing | 1.2.840.10008.5.1.4.1.1.1.1.1 | all from Table 23 | SCP | None |
| Enhanced MR Image Storage | 1.2.840.10008.5.1.4.1.1.4.1 | all from Table 23 | SCP | None |
| Enhanced SR Storage | 1.2.840.10008.5.1.4.1.1.88.22 | all from Table 23 | SCP | None |
| Grayscale Softcopy Presentation State Storage | 1.2.840.10008.5.1.4.1.1.11.1 | all from Table 23 | SCP | None |
| Hardcopy Color Image Storage | 1.2.840.10008.5.1.1.30 | all from Table 23 | SCP | None |
| Hardcopy Grayscale Image Storage | 1.2.840.10008.5.1.1.29 | all from Table 23 | SCP | None |
| Key Object Selection Document | 1.2.840.10008.5.1.4.1.1.88.59 | all from Table 23 | SCP | None |
| MR Image Storage | 1.2.840.10008.5.1.4.1.1.4 | all from Table 23 | SCP | None |
| MR Spectroscopy Storage | 1.2.840.10008.5.1.4.1.1.4.2 | all from Table 23 | SCP | None |
| PET Image Storage | 1.2.840.10008.5.1.4.1.1.128 | all from Table 23 | SCP | None |
| Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.20 | all from Table 23 | SCP | None |
| RETIRED Nuclear Medicine Image Storage | 1.2.840.10008.5.1.4.1.1.5 | all from Table 23 | SCP | None |
| RETIRED Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6 | all from Table 23 | SCP | None |
| RETIRED Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3 | all from Table 23 | SCP | None |
| Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7 | all from Table 23 | SCP | None |
| Stored Print Storage | 1.2.840.10008.5.1.1.27 | all from Table 23 | SCP | None |
| Ultrasound Image Storage | 1.2.840.10008.5.1.4.1.1.6.1 | all from Table 23 | SCP | None |
| Ultrasound Multi-frame Image Storage | 1.2.840.10008.5.1.4.1.1.3.1 | all from Table 23 | SCP | None |
| X-ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 | all from Table 23 | SCP | None |
| X-ray Radiofluoroscopic Image Storage | 1.2.840.10008.5.1.4.1.1.12.2 | all from Table 23 | SCP | None |

3.1.3.2.3 SOP Specific Conformance – Storage as SCP

DIAM⁴ Server conforms to the DICOM Storage Service Class at Level 2 (Full). No elements are discarded or coerced by DIAM⁴ Server. DIAM⁴ Server returns one of the following status codes (Table 34).

Table 34 : Storage status codes

| Service Status | Further Meaning | Protocol Codes | Related Fields | Description |
|----------------|-----------------------------------|----------------|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Refused | Out of resources | A700 | | Indicates that there was not enough storage space to store the image. Recovery from this condition is left to the administrative functions available in. |
| | SOP Class not supported | A800 | | Indicates that the SOP Class of the Image in the C-Store operation did not match the Abstract Syntax negotiated for the Presentation Context. |
| Error | Data set does not match SOP Class | A900 | | Indicates that the Data Set does not encode an instance of the SOP Class specified. |
| | Failed | C000 | | The operation was not successful. |
| | Cannot understand | C005 | | Indicates that the Data Set cannot be parsed into elements by |
| Warning | Coercion of data elements | B000 | | Data elements were modified before being stored. |
| | Data set does not match SOP Class | B007 | | Indicates that the Data Set does not match the SOP Class, but that the image was stored anyway. |
| | Elements Discarded | B006 | | Indicates that some of the elements of the Data Set were discarded. |
| | Duplicate SOP Instance UID | D000 | | Indicates that the SOP Instance UID of the specified image is already stored in the database. |
| Success | Success | 0000 | | Operation performed properly. |

3.1.3.2.4 Presentation Context Acceptance Criterion – Storage as SCP

DIAM⁴ Server will accept any number of Storage Presentation Contexts per association request. Any Abstract Syntax may be specified more than once in an association request, if the Transfer Syntaxes differ between the Presentation Contexts.

3.1.3.3 Real World Activity – Find

3.1.3.3.1 Associated Real World Activity – Find

DIAM⁴ Server will respond to query requests that are sent to it from an SCU.

3.1.3.3.2 Presentation Context Table – Find

DIAM⁴ Server will accept any of the transfer syntaxes listed in Table 35. For a QUERY request, DIAM⁴ Server will accept any of the Presentation Contexts listed in Table 36.

Table 35 : Find Transfer Syntaxes

| Transfer Syntax | UID |
|-----------------|-----|
|-----------------|-----|

| | |
|-------------------------------------------------|---------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |
| DICOM Explicit VR Big Endian Transfer Syntax | 1.2.840.10008.1.2.2 |

Table 36 : Find Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|-------------------------------------------|-----------------------------|-------------------|------|----------------------|
| SOP Class | SOP Class UID | | | |
| Patient Root Query/Retrieve IM Find | 1.2.840.10008.5.1.4.1.2.1.1 | all from Table 35 | SCP | See Note 1 |
| Study Root Query/Retrieve IM Find | 1.2.840.10008.5.1.4.1.2.2.1 | all from Table 35 | SCP | See Note 1 |
| Patient Study Only Query/Retrieve IM Find | 1.2.840.10008.5.1.4.1.2.3.1 | all from Table 35 | SCP | See Note 1 |

Note 1: Find Extended Negotiation will be supported. DIAM⁴ Server will negotiate with the following information:

Table 37 : Find Extended Negotiation

| Field Name | Value | Description of Field |
|--------------------|-------|------------------------------|
| Relational queries | 1 | Relational queries supported |

3.1.3.3.3 SOP Specific Conformance – Find

SOP classes of the Query/Retrieve Service Class are implemented via the DIMSE C-FIND and C-MOVE services as defined in Part 7 of the DICOM standard.

DIAM⁴ Server supports hierarchical queries. DIAM⁴ Server supports relational queries. DIAM⁴ Server, by default, supports all mandatory search keys.

Table 38 : Patient level attributes

| Description | Tag |
|--------------|------------------|
| Patient name | (0x0010, 0x0010) |
| Patient id | (0x0010, 0x0020) |

Table 39 : Study level attributes

| Description | Tag |
|----------------------------|------------------|
| Study id | (0x0020, 0x0010) |
| Study date | (0x0008, 0x0020) |
| Study time | (0x0008, 0x0010) |
| Accession number | (0x0008, 0x0050) |
| Modalities in Study | (0x0008, 0x0061) |
| Referring Physician's Name | (0x0008, 0x0090) |
| Study Description | (0x0008, 0x1030) |
| Study Instance UID | (0x0020, 0x000D) |

Table 40 : Series level attributes

| Description | Tag |
|---------------------|------------------|
| Series instance UID | (0x0020, 0x000E) |
| Series number | (0x0020, 0x0011) |
| Modality | (0x0008, 0x0060) |
| Series Description | (0x0008, 0x103e) |

| | |
|-------------------------------------|------------------|
| Request Attribute Sequence | (0x0040, 0x0275) |
| > Requested procedure ID | (0x0040, 0x1001) |
| > Scheduled Procedure Step ID | (0x0040, 0x0009) |
| Performed Procedure Step Start Date | (0x0040, 0x0244) |
| Performed Procedure Step Start Time | (0x0040, 0x0245) |

Table 41 : Composite Object Instance Level

| Description | Tag |
|------------------|------------------|
| SOP instance UID | (0x0008, 0x0018) |
| SOP Class UID | (0x0008, 0x0016) |
| Instance Number | (0x0020, 0x0013) |

Table 42 : SR Document Specific Level

| Description | Tag |
|-----------------------------|------------------|
| Completion Flag | (0x0040, 0xA491) |
| Verification Flag | (0x0040, 0xA493) |
| Verifying Observer Sequence | (0x0040, 0xA073) |
| > Verification DateTime | (0x0040, 0xA030) |
| > Verifying Observer Name | (0x0040, 0xA075) |
| Concept Name Code Sequence | (0x0040, 0xA043) |
| > Code Value | (0x0008, 0x0100) |
| > Coding Scheme Designator | (0x0008, 0x0102) |

Table 43 : Key Object Document Specific Level

| Description | Tag |
|----------------------------|------------------|
| Concept Name Code Sequence | (0x0040, 0xA043) |
| > Code Value | (0x0008, 0x0100) |
| > Coding Scheme Designator | (0x0008, 0x0102) |

Table 44 : Presentation State Specific Level attributes

| Description | Tag |
|-------------------------------|------------------|
| Content Label | (0x0070, 0x0080) |
| Content Description | (0x0070, 0x0081) |
| Presentation Creation Date | (0x0070, 0x0082) |
| Presentation Creation time | (0x0070, 0x0083) |
| Presentation Creator's Name | (0x0070, 0x0084) |
| Referenced Series Sequence | (0x0008, 0x1115) |
| > Series Instance UID | (0x0020, 0x000E) |
| > Referenced Image Sequence | (0x0008, 0x1140) |
| >>Referenced SOP Class UID | (0x0008, 0x1150) |
| >>Referenced SOP Instance UID | (0x0008, 0x1155) |

DIAM⁴ Server returns one of the following status codes to a C-FIND request.

Table 45 : C-FIND status codes

| Service Status | Further Meaning | Protocol Codes | Related Fields | Description |
|-----------------------|-------------------------------------------|-----------------------|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Refused | Out of Resources | A700 | | |
| Failed | Identifier does not match SOP Class | A900 | | The specified identifier contains a request that does not match the specified SOP Class. |
| | Unable to process | C001 | | For some reason (database off-line?) we cannot process this request at this time. |
| Cancel | Matching terminated due to Cancel Request | FE00 | | The original requester canceled this operation. |
| Pending | Pending | FF00 | | All Optional Keys are supported in the same manner as Required Keys. |
| | Pending | FF01 | | The matching operation is continuing. Warning that one or more Optional Keys were not supported in the same manner as Required Keys. |
| Success | Success | 0000 | | Operation performed properly. |

3.1.3.3.4 Presentation Context Acceptance Criterion - Find

DIAM⁴ Server will accept any number of Find Presentation Contexts per association request. Any Abstract Syntax may be specified more than once in an association request, if the Transfer Syntaxes differ between the Presentation Contexts.

3.1.3.4 Real World Activity – Move

3.1.3.4.1 Associated Real World Activity - Move

DIAM⁴ Server will respond to retrieve requests that are sent to it from an SCU.

3.1.3.4.2 Presentation Context Table – Move

DIAM⁴ Server will accept any of the transfer syntaxes listed in Table 46. For a MOVE request, DIAM⁴ Server will accept any of the Presentation Contexts listed in Table 47.

Table 46 : Move Transfer Syntaxes

| Transfer Syntax | UID |
|-------------------------------------------------|---------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |
| DICOM Explicit VR Big Endian Transfer Syntax | 1.2.840.10008.1.2.2 |

Table 47 : Move Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|-------------------------------------------|-----------------------------|-------------------|------|----------------------|
| SOP Class | SOP Class UID | | | |
| Patient Root Query/Retrieve IM Move | 1.2.840.10008.5.1.4.1.2.1.2 | all from Table 45 | SCP | See note 2 |
| Study Root Query/Retrieve IM Move | 1.2.840.10008.5.1.4.1.2.2.2 | all from Table 45 | SCP | See note 2 |
| Patient/Study Only Query/Retrieve IM Move | 1.2.840.10008.5.1.4.1.2.3.2 | all from Table 45 | SCP | See note 2 |

Note 2: Move Extended Negotiation will be supported. DIAM⁴ Server will negotiate with the following information:

Table 48 : Move Extended Negotiation

| Field Name | Value | Description of Field |
|----------------------|-------|--------------------------------|
| Relational-retrieval | 1 | Relational retrieval supported |

3.1.3.4.3 SOP Specific Conformance - Move

DIAM⁴ Server will try to establish an association with the move destination specified in the Move request. One or more of the Presentation Contexts listed in the Store section of this document may be negotiated in this association.

DIAM⁴ Server returns one of the following status codes to a C-MOVE request.

Table 49 : C-MOVE status codes

| Service Status | Further Meaning | Protocol Codes | Related Fields | Description |
|----------------|------------------------------------------|----------------|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Refused | Out of Resources | A701 | | Unable to calculate number of matches. |
| | Out of Resources | A702 | | Unable to perform storage of images to move destination. |
| Failed | Move destination unknown | A801 | | The destination of this move request is unknown. |
| | Identifier does not match SOP Class | A900 | | The specified identifier contains a request that does not match the specified SOP Class. |
| | Unable to process | C002 | | Indicates that Diam4 cannot process this request at this time. |
| Cancel | Storage terminated due to Cancel Request | FE00 | | The original requester canceled this operation. |
| Warning | Warning | B000 | | Storage complete with one or more failures. |
| Pending | Pending | FF00 | | The storage operation is continuing. |
| | Pending for a long time | FF02 | | This operation is expected to require a long period of time to complete. The SCU may break the association at any time, but the operation will continue to completion. |
| Success | Success | 0000 | | Operation performed properly. |

3.1.3.4.4 Presentation Context Acceptance Criterion - Move

DIAM⁴ Server will accept any number of Move Presentation Contexts per association request. Any Abstract Syntax may be specified more than once in an association request, if the Transfer Syntaxes differ between the Presentation Contexts.

3.1.3.4.5 Transfer Syntax Selection Policies - Move

By default, DIAM⁴ Server sends the IOD using the transfer syntax that was used when the image was originally stored.

3.1.3.5 Real World Activity – Storage Commitment

3.1.3.5.1 Associated Real World Activity – Storage Commitment

DIAM⁴ Server stores images that are sent to it from an *SCU*. The request for storage commitment may then be transmitted to DIAM⁴ Server together with a list of references to one or more SOP instances. DIAM⁴ Server will receive and respond to DIMSE N-ACTION. The following message is supported:

- Request Storage Commitment - to request the safekeeping of a set of SOP instances

3.1.3.5.2 Presentation Context Table – Storage Commitment

DIAM⁴ Server will accept any of the transfer syntaxes listed in Table 50. For a Storage commitment request, DIAM⁴ Server will accept any of the Presentation Contexts listed in Table 51.

Table 50 : Transfer Syntaxes

| Transfer Syntax | UID |
|-------------------------------------------------|---------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |
| DICOM Explicit VR Big Endian Transfer Syntax | 1.2.840.10008.1.2.2 |

Table 51 : Storage Commitment Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|---------------------------------------------------|----------------------------------------------|-------------------|------|----------------------|
| SOP Class Storage Commitment Push Model | SOP Class UID 1.2.840.10008.1.20.1 | all from Table 50 | SCP | None |

3.1.3.5.3 SOP Specific Conformance – Storage Commitment

DIAM⁴ Server supports the following elements for this SOP class as an SCP:

Table 52: Storage Commitment Request – Action Information

| Action Type Name | Action Type ID | Attribute Name | Tag |
|----------------------------|----------------|-------------------------------|-------------|
| Request Storage Commitment | 1 | Transaction UID | (0008,1195) |
| | | Referenced SOP Sequence | (0008,1199) |
| | | > Referenced SOP Class UID | (0008,1150) |
| | | > Referenced SOP Instance UID | (0008,1155) |

DIAM⁴ Server returns one of the following status codes:

Table 53 : Storage Commitment status codes

| Service Status | Further Meaning | Protocol Codes | Related Fields | Description |
|----------------|-----------------|----------------|----------------|-----------------------------------|
| Error | Failed | C000 | | The operation was not successful. |
| Success | Success | 0000 | | Operation performed properly. |

3.1.3.5.4 Result – Storage Commitment

If DIAM⁴ Server determines that it has successfully completed storage commitment, DIAM⁴ Server issues an N-EVENT-REPORT to the SCU including references to the successfully stored SOP Instances contained in the N-ACTION.

The N-EVENT-REPORT contains the Transaction UID value contained in the initiating N-ACTION. The N-EVENT-REPORT is sent on a separate association from the N-ACTION operation.

DIAM⁴ Server supports the Event Information as specified in Table 54.

Table 54: Storage Commitment Result – Event Information

| Action Type Name | Event Type ID | Attribute Name | Tag |
|---------------------------------------|---------------|------------------------------|-------------|
| Storage Commitment Request Successful | 1 | Transaction UID | (0008,1195) |
| | | Referenced SOP Sequence | (0008,1199) |
| | | >Referenced SOP Class UID | (0008,1150) |
| | | >Referenced SOP Instance UID | (0008,1155) |
| | | Failed SOP Sequence | (0008,1198) |
| | | >Referenced SOP Class UID | (0008,1150) |
| | | >Referenced SOP Instance UID | (0008,1155) |
| | | >Failure reason | (0008,1197) |

In a cache-only configuration, DIAM⁴ Server commits to storing a SOP Instance as long as there is available disk space. In this configuration, DIAM⁴ Server may delete SOP Instances.

3.1.3.6 Real World Activity – MPPS as SCP

3.1.3.6.1 Associated Real World Activity – Modality Performed

Procedure Step

If the functionality is activated in DIAM⁴ settings, the DIAM⁴ Server accepts MPPS messages from a remote AE. The messages are either N-CREATE or N-SET. The DIAM⁴ Server is able to treat many incoming messages in the same time.

When the DIAM⁴ Server receives a Modality Performed Procedure Step SOP Class N-CREATE request, it stores the information. When the DIAM⁴ Server receives a Modality Performed Procedure Step SOP Class N-SET request, it updates the information.

Furthermore, these messages are forwarded to the configured peer applications.

3.1.3.6.2 Presentation Context Table – Modality Performed

Procedure Step

DIAM⁴ Server will accept any of the transfer syntaxes listed in Table 55. DIAM⁴ Server will accept any of the Presentation Contexts listed in Table 56.

Table 55 : Transfer Syntaxes

| Transfer Syntax | UID |
|-------------------------------------------------|---------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |
| DICOM Explicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2.1 |

Table 56 : MPPS Presentation Contexts

| Abstract Syntax | | Transfer Syntax | Role | Extended Negotiation |
|-----------------------------------|-------------------------|------------------------|-------------|-----------------------------|
| SOP Class | SOP Class UID | | | |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | all from Table 55 | SCP | None |

3.1.3.6.3 SOP Specific Conformance – Modality Performed

Procedure Step

When receiving a MPPS N-CREATE request, DIAM⁴ Server creates a MPPS object with the "IN PROGRESS" status.

The content of the MPPS object is updated when receiving a MPPS N-SET request. The status is updated to either "IN PROGRESS" or "COMPLETED" or "DISCONTINUED".

If DIAM⁴ Server receives another MPPS N-SET request for this object it will accept if the current status in the database is "IN PROGRESS", else if this status is "COMPLETED" or "DISCONTINUED", it will return an error code as described in table 59.

DIAM4 Server supports the following elements MPPS N-CREATE request as an SCP :

Table 57 : MPPS N-CREATE data elements

| Description | Tag |
|-------------------------------------------------|-------------|
| Performed Procedure Step Relationship | |
| Patient's Name | (0010,0010) |
| Patient ID | (0010,0020) |
| Patient's Birth Date | (0010,0030) |
| Patient's Sex | (0010,0040) |
| Referenced Patient Sequence | (0008,1120) |
| > Referenced SOP Class UID | (0008,1150) |
| > Referenced SOP Instance UID | (0008,1155) |
| Scheduled Step Attributes Sequence | (0040,0270) |
| > Study Instance UID | (0020,000D) |
| > Referenced Study Sequence | (0008,1110) |
| >> Referenced SOP Class UID | (0008,1150) |
| >> Referenced SOP Instance UID | (0008,1155) |
| > Accession Number | (0008,0050) |
| > Placer Order Number / Imaging Service Request | (0040,2016) |
| > Filler Order Number / Imaging Service Request | (0040,2017) |
| > Requested Procedure ID | (0040,1001) |
| > Requested Procedure Description | (0032,1060) |
| > Scheduled Procedure Step ID | (0040,0009) |
| > Scheduled Procedure Step Description | (0040,0007) |
| > Scheduled Protocol Code Sequence | (0040,0008) |
| >> Code Value | (0008,0100) |
| >> Coding Scheme Designator | (0008,0102) |
| >> Coding Scheme Version | (0008,0103) |
| >> Code Meaning | (0008,0104) |
| Performed Procedure Step Information | |
| Performed Station AE Title | (0040,0241) |
| Performed Station Name | (0040,0242) |
| Performed Location | (0040,0243) |
| Performed Procedure Step Start Date | (0040,0244) |
| Performed Procedure Step Start Time | (0040,0245) |
| Performed Procedure Step ID | (0040,0253) |
| Performed Procedure Step End Date | (0040,0250) |
| Performed Procedure Step End Time | (0040,0251) |
| Performed Procedure Step Status | (0040,0252) |
| Performed Procedure Step Description | (0040,0254) |
| Comments on the Performed Procedure Step | (0040,0280) |
| Performed Procedure Type Description | (0040,0255) |

| | |
|---------------------------------------------------------------|-------------|
| Procedure Code Sequence | (0008,1032) |
| > Code Value | (0008,0100) |
| > Coding Scheme Designator | (0008,0102) |
| > Coding Scheme Version | (0008,0103) |
| > Code Meaning | (0008,0104) |
| Performed Procedure Step Discontinuation Reason Code Sequence | (0040,0281) |
| > Code Value | (0008,0100) |
| > Coding Scheme Designator | (0008,0102) |
| > Coding Scheme Version | (0008,0103) |
| > Code Meaning | (0008,0104) |
| Image Acquisition Results | |
| Modality | (0008,0060) |
| Study ID | (0020,0010) |
| Performed Protocol Code Sequence | (0040,0260) |
| > Code Value | (0008,0100) |
| > Coding Scheme Designator | (0008,0102) |
| > Coding Scheme Version | (0008,0103) |
| > Code Meaning | (0008,0104) |
| Performed Series Sequence | (0040,0340) |
| > Performing Physician's Name | (0008,1050) |
| > Operator's Name | (0008,1070) |
| > Protocol Name | (0018,1030) |
| > Series Instance UID | (0020,000E) |
| > Series Description | (0008,103E) |
| > Retrieve AE Title | (0008,0054) |
| > Referenced Image Sequence | (0008,1140) |
| >> Referenced SOP Class UID | (0008,1150) |
| >> Referenced SOP Instance UID | (0008,1155) |
| > Referenced Non-Image Composite SOP Instance Sequence | (0040,0220) |
| >> Referenced SOP Class UID | (0008,1150) |
| >> Referenced SOP Instance UID | (0008,1155) |
| All other attributes from Radiation Dose Module | |
| All other attributes from Billing and Material Code Module | |

DIAM⁴ Server supports the following elements MPPS N-SET request as an SCP:

Table 58 : MPPS N-SET data elements

| Description | Tag |
|---------------------------------------------------------------|-------------|
| Performed Procedure Step Information | |
| Performed Procedure Step End Date | (0040,0250) |
| Performed Procedure Step End Time | (0040,0251) |
| Performed Procedure Step Status | (0040,0252) |
| Performed Procedure Step Description | (0040,0254) |
| Comments on the Performed Procedure Step | (0040,0280) |
| Performed Procedure Type Description | (0040,0255) |
| Procedure Code Sequence | (0008,1032) |
| > Code Value | (0008,0100) |
| > Coding Scheme Designator | (0008,0102) |
| > Coding Scheme Version | (0008,0103) |
| > Code Meaning | (0008,0104) |
| Performed Procedure Step Discontinuation Reason Code Sequence | (0040,0281) |
| > Code Value | (0008,0100) |
| > Coding Scheme Designator | (0008,0102) |
| > Coding Scheme Version | (0008,0103) |
| > Code Meaning | (0008,0104) |
| Image Acquisition Results | |
| Performed Protocol Code Sequence | (0040,0260) |
| > Code Value | (0008,0100) |

| | |
|------------------------------------------------------------|-------------|
| > Coding Scheme Designator | (0008,0102) |
| > Coding Scheme Version | (0008,0103) |
| > Code Meaning | (0008,0104) |
| Performed Series Sequence | (0040,0340) |
| > Performing Physician's Name | (0008,1050) |
| > Operator's Name | (0008,1070) |
| > Protocol Name | (0018,1030) |
| > Series Instance UID | (0020,000E) |
| > Series Description | (0008,103E) |
| > Retrieve AE Title | (0008,0054) |
| > Referenced Image Sequence | (0008,1140) |
| >> Referenced SOP Class UID | (0008,1150) |
| >> Referenced SOP Instance UID | (0008,1155) |
| > Referenced Non-Image Composite SOP Instance Sequence | (0040,0220) |
| >> Referenced SOP Class UID | (0008,1150) |
| >> Referenced SOP Instance UID | (0008,1155) |
| All other attributes from Radiation Dose Module | |
| All other attributes from Billing and Material Code Module | |

DIAM⁴ Server returns one of the following status codes.

Table 59 : MPPS N-SET status codes

| Service Status | Further Meaning | Protocol Codes | Related Fields | Description |
|-----------------------|-------------------------------------------------------------------------------|-----------------------|-----------------------|---------------------------------------------------------------------------|
| Failure | Processing Failure – Performed Procedure Step Object may no longer be updated | 0110 | | No update can be performed if the status is "COMPLETED" or "DISCONTINUED" |
| Success | Success | 0000 | | Operation performed properly. |

3.1.4 Presentation State Display

DIAM⁴ Server manages presentation states. It applies them to the related images when these images are displayed. The following specific presentation modules are supported:

Table 60 : Presentation State Module (C.11.10)

| Attribute Name | Tag | Comments |
|-------------------------------|-------------|----------|
| Instance Number | (0020,0013) | |
| Content Label | (0070,0080) | |
| Content Description | (0070,0081) | |
| Presentation Creation Time | (0070,0083) | |
| Content Creator's Name | (0070,0084) | |
| >Series Instance UID | (0020,000E) | |
| >Referenced Image Sequence | (0008,1140) | |
| >>Referenced SOP Class UID | (0008,1150) | |
| >>Referenced SOP Instance UID | (0008,1155) | |
| >>Referenced Frame Number | (0008,1160) | |
| Shutter Presentation Value | (0018,1622) | |

Table 61 : Display Shutter Module (C.7.6.11)

| Attribute Name | Tag | Comments |
|-----------------------------------|-------------|----------|
| Shutter Shape | (0018,1600) | |
| Shutter Left Vertical Edge | (0018,1602) | |
| Shutter Right Vertical Edge | (0018,1604) | |
| Shutter Upper Horizontal Edge | (0018,1606) | |
| Shutter Lower Horizontal Edge | (0018,1608) | |
| Center of Circular Shutter | (0018,1610) | |
| Radius of Circular Shutter | (0018,1612) | |
| Vertices of the Polygonal Shutter | (0018,1620) | |
| Shutter Presentation Value | (0018,1622) | |

N.B.: the Bitmap Display Shutter Module is not supported

Table 62 : Overlay Plane (C.9.2)

| Attribute Name | Tag | Comments |
|------------------------|-------------|------------------------------------------------|
| Overlay Rows | (60xx,0010) | |
| Overlay Columns | (60xx,0011) | |
| Overlay Type | (60xx,0040) | Enumerated Values: G = Graphics R = ROI. |
| Overlay Origin | (60xx,0050) | |
| Overlay Bits Allocated | (60xx,0100) | |
| Overlay Bit Position | (60xx,0102) | |
| Overlay Data | (60xx,3000) | |
| Overlay Description | (60xx,0022) | |
| Overlay Label | (60xx,1500) | |

Table 63 : Overlay/Curve Activation (C.11.7)

The activation of overlays is not supported: DIAM4 Server displays either all overlays of the image referenced, or none of them.

Table 64 : Displayed Area (C.10.4)

| Attribute Name | Tag | Comments |
|------------------------------------------|-------------|------------------------------------------------------------|
| Displayed Area Selection Sequence | (0070,005A) | |
| >Referenced Image Sequence | (0008,1140) | |
| >>Referenced SOP Class UID | (0008,1150) | |
| >>Referenced SOP Instance UID | (0008,1155) | |
| >>Referenced Frame Number | (0008,1160) | |
| >Displayed Area Top Left Hand Corner | (0070,0052) | |
| >Displayed Area Bottom Right Hand Corner | (0070,0053) | |
| >Presentation Size Mode | (0070,0100) | Enumerated Values: SCALE TO FIT TRUE SIZE MAGNIFY |
| >Presentation Pixel Spacing | (0070,0101) | |
| >Presentation Pixel Aspect Ratio | (0070,0102) | |
| >Presentation Pixel Magnification Ratio | (0070,0103) | |

Table 65 : Graphic Annotation (C.10.5)

| Attribute Name | Tag | Comments |
|----------------------------------------------|-------------|-------------------------------------|
| Graphic Annotation Sequence | (0070,0001) | |
| >Referenced Image Sequence | (0008,1140) | |
| >>Referenced SOP Class UID | (0008,1150) | |
| >>Referenced SOP Instance UID | (0008,1155) | |
| >>Referenced Frame Number | (0008,1160) | |
| >Graphic Layer | (0070,0002) | |
| >Text Object Sequence | (0070,0008) | |
| >>Bounding Box Annotation Units | (0070,0003) | |
| >>Anchor Point Annotation Units | (0070,0004) | |
| >>Unformatted Text Value | (0070,0006) | |
| >>Bounding Box Top Left Hand Corner | (0070,0010) | |
| >>Bounding Box Bottom Right Hand Corner | (0070,0011) | |
| >>Bounding Box Text Horizontal Justification | (0070,0012) | |
| >>Anchor Point | (0070,0014) | |
| >>Anchor Point Visibility (0070,0015) | (0070,0015) | |
| >Graphic Object Sequence | (0070,0009) | |
| >>Graphic Annotation Units | (0070,0005) | |
| >>Graphic Dimensions | (0070,0020) | |
| >>Number of Graphic Points | (0070,0021) | |
| >> Graphic Data | (0070,0022) | |
| >>Graphic Type | (0070,0023) | Type INTERPOLATED is not supported. |
| >>Graphic Filled | (0070,0024) | |

Table 66 : Spatial Transformation (C.10.6)

| Attribute Name | Tag | Comments |
|-----------------------|-------------|----------|
| Image Rotation | (0070,0041) | |
| Image Horizontal Flip | (0070,0041) | |

Table 67 : Modality LUT (C.11.1)

| Attribute Name | Tag | Comments |
|-----------------------|-------------|----------|
| Modality LUT Sequence | (0028,3000) | |
| >LUT Descriptor | (0028,3002) | |
| >LUT Explanation | (0028,3003) | |
| >Modality LUT Type | (0028,3004) | |
| >LUT Data | (0028,3006) | |
| Rescale Intercept | (0028,1052) | |
| Rescale Slope | (0028,1053) | |
| Rescale Type | (0028,1054) | |

Table 68 : Softcopy VOI LUT (C.11.8)

| Attribute Name | Tag | Comments |
|------------------------------------|-------------|----------|
| Softcopy VOI LUT Sequence | (0028,3110) | |
| >Referenced Image Sequence | (0008,1140) | |
| >>Referenced SOP Class UID | (0008,1150) | |
| >>Referenced SOP Instance UID | (0008,1155) | |
| >>Referenced Frame Number | (0008,1160) | |
| >VOI LUT Sequence | (0028,3010) | |
| >>LUT Descriptor | (0028,3002) | |
| >>LUT Explanation | (0028,3003) | |
| >>LUT Data | (0028,3006) | |
| >Window Center | (0028,1050) | |
| >Window Width | (0028,1051) | |
| >Window Center & Width Explanation | (0028,1055) | |

Table 69 : Softcopy Presentation LUT (C.11.6)

| Attribute Name | Tag | Comments |
|---------------------------|-------------|-----------------|
| Presentation LUT Sequence | (2050,0010) | |
| >LUT Descriptor | (0028,3002) | |
| >LUT Explanation | (0028,3003) | |
| >LUT Data | (0028,3006) | |
| Presentation LUT Shape | (2050,0020) | |

4 Communications Profiles

DIAM⁴ Server provides DICOM V3.0 TCP/IP Network Communication Support as defined in Part 8 of the DICOM Standard.

4.1 TCP/IP Stack

DIAM⁴ Server inherits its TCP/IP stack from the computer system upon which it executes.

4.1.1 Physical Media Support

DIAM⁴ Server is indifferent to the physical medium over which TCP/IP executes; it inherits the medium from the computer system upon which it executes.

5 Extensions/Specializations/Privatizations

6 Configuration

DIAM⁴ Server obtains configuration information from the following sources:

Mapping from Application Entity Title to Presentation Address is provided by the database. Along with this mapping, the database stores those AE titles that are allowed to communicate with DIAM⁴ Server.

7 Support for Extended Character Sets

DIAM⁴ Server is known to support the following character sets:

| | |
|--------------------|----------------------|
| ISO-IR 6 (default) | Basic G0 Set |
| ISO-IR 100 | Latin Alphabet No. 1 |