



# Global Imaging Online

WEB PACS, WEB RIS, Télé-Imagerie, Applications 3D

## IHE Integration Statement 18 Mai 2010

Vendor Global Imaging Online  
Product name DIAM 4  
Release 3.0

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Image Manager / Archive	None
	Report Reader	None
	Report Repository	None
Audit Trail and Node Authentication	Secure Application	None
Cardiac Catheterization Workflow	Image Manager/Archive	Intermittently Connected Modality
	Performed Procedure Step Manager	None
Consistent Presentation of Images	Image Manager/Archive	None
Consistent Time	Time Client	None
Cross-Enterprise Clinical Documents Share (XDS.a)	Imaging Document Source	None
Cross-enterprise Document Sharing for Imaging	Document Source	None
Echocardiography Workflow	Image Manager/Archive	Intermittently Connected Modality
	Performed Procedure Step Manager	None

Profiles Implemented	Actors Implemented	Options Implemented
Evidence Documents	Image Manager/Archive	None
Key Image Note	Image Manager/Archive	None
Mammography Acquisition Workflow	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Mammography Image	Image Manager/Archive	None
Patient Information Reconciliation	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Portable Data for Imaging	Image Display	None
	Portable Media Creator	Web Content, UMedia for DPI, DVD Media
	Portable Media Importer	None
Scheduled Workflow	Image Manager/Archive	None
	Performed Procedure Step Manager	None
Simple Image and Numeric Report	Report Reader	None
	Report Repository	None

<b>Internet address for vendor's IHE information</b> <a href="http://www.global-imaging.net/conformance.htm">http://www.global-imaging.net/conformance.htm</a>	
Links to Standards Conformance Statements for the implementation	
HL7	<a href="http://www.global-imaging.net/conformance.htm">http://www.global-imaging.net/conformance.htm</a>
DICOM:	

Links to general information on IHE

In North America : <a href="http://www.ihe.net">http://www.ihe.net</a>	<b>In Europe :</b> <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	<b>In Japan :</b> <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>
---	---	--