



Global Imaging Online

WEB PACS, WEB RIS, Télé-Imagerie, Applications 3D

IHE Integration Statement 20 Mai 2010

Vendor Global Imaging Online
Product name: DIAM Box EFI
Release 1.6

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Consistent Presentation of Images	Image Display	None
	Print Composer	None
Evidence Documents	Evidence Creator	None
	Image Display	None
Key Image Note	Evidence Creator	None
	Image Display	None
Mammography Image	Image Display	None
	Print Composer	None
Scheduled Workflow	Evidence Creator	None
	Image Display	None

Internet address for vendor's IHE information http://www.global-imaging.net/conformance.htm	
Links to Standards Conformance Statements for the implementation	
HL7 DICOM	http://www.global-imaging.net/conformance.htm

Links to general information on IHE

In North America :http://www.ihe.net	In Europe : http://www.ihe-europe.net	In Japan : http://www.ihe-j.org
--	---	--